

Unit-02 Basic Measurement II

Objective :

Use spherometer to measure the curvature radius of a spherical surface.

Apparatus :

Spherometer, traditional spherometer, plate glass, convex glass

Principle :

A. The Construction

The structure of a spherometer is shown as Fig.1. **A** is a tripod with three arms **B**, **C**, and **D**, and that also with equivalent length. The included angle of each two arms is 120° . **E** is the main-meter, one tick is 1 mm. And there are 100 ticks on the periphery of sub-meter **F**. Probe **H** shifts 1 mm when sub-meter is screwed one round, which means the length per tick is 0.01 mm. We can get the height of probe from the main-meter and sub-meter.




Figure 1. The structure of a spherometer

B. Measurement

(a) Main-ruler reading

Sub-ruler's pointer between N and $N+1$ mm, the readings record N mm.

(b) Sub-ruler reading

Find out the ticks on sub-meter. It should include estimate value.

C. Curvature Radius

As shown in Fig.2, the distance between each arm of the tripod is S , hence **B**, **C** and **D** form an equilateral triangle. We can make a circumscribed circle of radius r .

The extended line of probe **H** is sure to pass the center of the sphere **O'**. We assume that the line **O'H** intersects plane **BCD** at **O**. Let curvature radius is R , if we measure the height is a by probe **H**, and then we could calculate curvature radius R by Pythagorean proposition.

$$r = \frac{S}{\sqrt{3}}$$

$$R^2 = (R - a)^2 + r^2$$

$$R = \frac{S^2}{6a} + \frac{a}{2}$$




Figure 2. Vertical view and Side view of a spherometer

Remarks :

1. When you put spherometer on the glossy glass plate, please gently lay down.

Procedure :

1. Adjust the fixed tips to make it into triangle.
2. Push down 3 tips B, C, D on the paper.
3. Remove the spherometer and draw a triangle.
4. Use vernier caliper to measure the length between two fixed tips.
5. Put the spherometer on the plate glass, and making the tips of **BCDH** contact the plate glass. Recording the value a_0 .
6. Put the spherometer on the convex glass, and making the tips of **BCDH** contact the convex glass. Recording the value a_1 .

7. The difference between a_0 and a_1 is a , that is the height of **H** from plane **BCD**.

$$a = |a_1 - a_0|$$

8. Get the average value and standard deviation of the mean.
9. Calculate the curvature radius R . (should consider the error transfer)

[Note] Curvature Radius $R = \bar{R} \pm \sigma_R$

$$\bar{R} = \frac{\bar{S}^2}{6a} + \frac{\bar{a}}{2}$$

$$\sigma_R = \sqrt{\sigma_s^2 \left(\frac{\bar{S}}{3a} \right)^2 + \sigma_a^2 \left(-\frac{\bar{S}^2}{6a^2} + \frac{1}{2} \right)^2}$$

10. Used traditional spherometer to do this experiment.

Questions :

1. The construction of traditional spherometer and spherometer are the same? Please explain.
2. Can we use spherometer to measure the radius of convex len? Please explain.